






## Compatibility of Supraglottic Airway Devices with ETTs

SDA\ETT	Sizes	6.0	6.5	7	7.5	8	8.5
<b>I-GEL</b> 	3	yes	no	no	no	no	NA
	4	yes	yes	yes	no	no	NA
	5	yes	yes	yes	yes	yes	NA
<b>Ambu</b> 	3	yes	yes	no	no	no	NA
	4	yes	yes	yes	yes	no	NA
	5	yes	yes	yes	yes	yes	NA
<b>Classic</b> 	3	yes	no	no	no	no	NA
	4	yes	no	no	no	no	NA
	5	yes	yes	yes	no	no	NA
<b>Air-Q</b> 	2.5	yes	yes	no	no	no	no
	3.5	yes	yes	yes	yes	no	no
	4.5	yes	yes	yes	yes	yes	yes
<b>Supreme</b> 	All sizes	<b>Non-Compatible with ETT</b>					

SAD: Supraglottic Airway Device

ETT: Endotracheal Tube

Appreciate to have your feedback sent to: [KACCS@kma.org.kw](mailto:KACCS@kma.org.kw)